English

Quinn's Advantage™ Protein Plus Fertilization (HTF) Medium

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Protein Plus Fertilization (HTF) Medium	ART-1520	20 mL

INTENDED USE

This product is intended for in vitro procedures involving the fertilization of human oocytes.

PRODUCT DESCRIPTION

This medium is a modification of Human Tubal Fluid (HTF) that was first described by Quinn and colleagues (Fertil Steril. 1984;41:202, 1985;44:493). Modifications to the original formulation include the addition of citrate, selected nonessential amino acids, taurine, and the presence of lactate in the form of calcium lactate.

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

To avoid problems with contamination, handle using aseptic techniques and discard any

excess product that remains in the bottle or vial after procedure is completed.

This product contains 3mg/mL human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation / removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Single use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) end-ing with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

This is the preferred medium for fertilization procedures in vitro. We currently recommend the use of Quinn's Advantage™ Protein Plus Cleavage Medium (REF # ART-1526) for culture of embryos from Day 1 until Day 3, and Quinn's Advantage™ Protein Plus Blastocyst Medium (REF # ART-1529) for culture from Day 3 to Day 5/6. The media require no protein supplement.

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature and equilibrate with desired atmosphere containing 5% CO₂ prior to use. Do not freeze or expose to temperatures greater than 39 °C. It is advisable to regass the medium with 5% CO₂ before replacing opened bottles back in the refrigerator.

In order to ensure optimal performance, we strongly recommend measuring the pH of the medium under laboratory working conditions and adjusting the level of CO₂ used to attain the desired pH range for optimal embryo development. The desired pH range for Quinn's Advantage™ Protein Plus Fertilization (HTF) Medium is 7.3 ± 0.1.

The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

ART-1526 Quinn's Advantage™ Protein Plus Cleavage Medium

ART-1529 Quinn's Advantage™ Protein Plus Blastocyst Medium

SAGE™ In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is a trademark of CooperSurgical, Inc.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

EXPLANATION OF SYMBOLS

- REF** Catalog Number
- LOT** Batch Number
- Use By (year, month, day)
- Do Not Reuse
- Temperature Limitation
- Aseptic Technique Sterilization Membrane Filtered (SAL 10⁻³)
- ATTENTION:** See instructions for use.
- EC REP** Authorized Representative in the European Community.
- Product conforms to the Medical Device Directive 93/42/EEC
- Manufacturer
- RX ONLY** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

SAGE In Vitro Fertilization
a CooperSurgical Company

SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

EC REP
ORIGIO a/s
Knardrupvej 2
DK-2760 Måløv
Denmark
www.origio.com

Customer Service:
E-mail: customer.service@origio.com
Tel: +45 46 79 02 02 Fax: +45 46 79 03 02

Français

Quinn's Advantage™ Protein Plus Fertilization (HTF) Medium

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Quinn's Advantage™ este o marcă comercială a CooperSurgical, Inc.