English

Quinn's™ Sperm Washing Medium

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's™ Sperm Washing Medium	ART-1005	12 x 12 mL
Quinn's™ Sperm Washing Medium	ART-1006	100 mL

INTENDED USE

This product was intended for in vitro procedures involving sperm washing and incubation (Seminars Reprod Endocrinol. 1987;5:23).

PRODUCT DESCRIPTION

This medium is a modification of HEPES-buffered Human Tubal Fluid (HTF) that was first described by Quinn and colleagues (Fertil Steril. 1984;41:202, 1985;44:493).

This product contains 10 mg/L of genta-micin, an aminoglycoside antibiotic.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

Quinn's™ Sperm Washing Medium should be tightly capped when used in a CO₂ incubator to avoid pH levels of 7.0 or less.

The Quinn's™ Sperm Washing Medium (ART-1005 and ART-1006) contains 5 mg/mL of human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Single use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if

the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic genta-micin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

Information on sperm processing procedures for intrauterine insemination (IUI) can be obtained by contacting the SAGE™ Support Line. Quinn's™ Sperm Washing Medium can be used with the PureCep™ line of sperm separation products.

BUFFER SYSTEM

Quinn's™ Sperm Washing Medium uses a buffering system composed of a 21 mM HEPES (N-2-Hydroxy-ethylpiperazine-N1-2-ethanesulfonic acid) and 4 mM sodium bicarbonate combination. This buffering system provides optimum pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO₂ incubator.

Each laboratory should make its own determination of which medium to use for each particular procedure.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to ambient or incubator (37 °C) temperature prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- A. Remove desired volume of product using aseptic procedures.
- B. Once removed, do not return any volume of product to the original container.
- C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

SAGE™ In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's™ is a trademark of CooperSurgical, Inc.

PureCep™ is a trademark of CooperSurgical, Inc.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

EXPLANATION OF SYMBOLS

-  Catalog Number
-  Batch Number
-  Use By (year, month, day)
-  Do Not Reuse
-  Temperature Limitation
-  Aseptic Technique Sterilization
Membrane Filtered (SAL 10⁻³)
-  **ATTENTION:**
See instructions for use.
-  Authorized Representative
in the European Community.
-  Product conforms to the
Medical Device Directive
93/42/EEC
-  Manufacturer
-  **RX ONLY**
U.S. Federal law restricts this
device to sale by or on the order
of a physician (or properly licensed
practitioner).

In Vitro Fertilization
a CooperSurgical Company

SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

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Quinn's™ este o marcă comercială a CooperSurgical, Inc.
PureCception™ este o marcă comercială a CooperSurgical, Inc.