5652-03
Ver. 3: 2016.Aug.24
1/1

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en - english

CSC (Choline Substituted Cryopreservation) Thawing Medium

CSC Thawing Medium Kit
For laboratory procedures only.

CSC Thawing Medium Kit

ART-8018

Sufficient for four (4) thawing procedures

INTENDED USE

This product is intended for thawing of mature metaphase II oocytes frozen with ART-8017, CSC Freezing Medium. This kit is designed to be used in conjunction with the SAGE™ CSC Freezing Medium Kit (ART-8017).

PRODUCT DESCRIPTION
CSC Thawing Medium will allow for the efficient thawing of mature metaphase II oocytes. The components and recommended procedures are the preferred method for improved oocyte survivability.

CONTAIN:
10 mg/L of gentamicin, an aminoglycoside antibiotic.
12 mg/mL human serum albumin.

MATERIALS PROVIDED IN THE THAWING KIT
1. 4 x 2 mL vials of CSC Thawing Medium 0.5 M Sucrose (REF # ART-8018-A)
2. 4 x 2 mL vials of CSC Thawing Medium 0.2 M Sucrose (REF # ART-8018-B)
3. 4 x 4 mL vials of CSC Washing Medium HEPES Buffered HTF (REF # ART-8018-C)

PRECAUTIONS AND WARNINGS
Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

To avoid problems with contamination, use aseptic techniques and discard any excess product that remains in the vial after procedure is completed.

Caution: All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HBsAg, and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore source material have been tested for parvovirus B19 and found non-elevated. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.
This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that that patient is not sensitized to this antibiotic.

Single use: To avoid problems with contamination, use aseptic techniques and discard any excess product that remains in the vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

QUALITY ASSURANCE
One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE FOR THAWING OF OOCYTES
1. Pre-procedure Set-up: Remove one vial each of CSC Thawing Medium 0.5 M Sucrose (REF # ART-8018-A), CSC Thawing Medium 0.2 M Sucrose (REF # ART-8018-B), and CSC Washing Medium HEPES Buffered HTF (REF # ART-8018-C) from storage at 2–8 °C and warm to room temperature (22-24 °C).
2. For each thawing procedure, set up four dishes. The first dish contains 2 mL of the CSC Thawing Medium 0.5 M Sucrose, the second dish contains 2 mL of the CSC Thawing Medium 0.2 M Sucrose, and the 3rd and 4th dish both contain 2 mL of the CSC Washing Medium. It is recommended that the media be covered with the preferred oil during use to minimize evaporation of water and a subsequent change in osmolality of the media.

THAWING PROTOCOL
1. If oocytes were frozen in straws, they should be thawed rapidly (at least 275 °C/min). Initially hold the straw in air for 30 to 40 seconds and then immerse it in a water bath at 30-35 °C until the ice has fully melted. Cryovials are thawed by immersion in a 32 °C water bath until all ice crystals have disappeared. Thaw only one cryovial at a time.
2. Transfer the liquid contents of the thawed straw or cryovial to a dry dish and quickly locate the oocytes. Pick the oocytes up in a minimal amount of medium and transfer them first to 2 mL of CSC Thawing Medium 0.5 M Sucrose (REF # ART-8018-A) at ambient temperature (22-24 °C) for 10 minutes.
3. Then, using a fine-bore pipette, transfer the oocytes from the CSC Thawing Medium 0.5 M Sucrose to 2 mL of CSC Thawing Medium 0.2 M Sucrose (REF # ART-8018-B) at ambient temperature for 10 minutes. Use a new transfer pipette for each procedure to minimize the carry-over of cryoprotectant from one solution to the next.
4. The oocytes are then washed through two dishes containing 2 mL of CSC Washing Medium (REF # ART-8018-C) at 37 °C for 5 minutes in each solution. The oocytes are thoroughly washed by pipetting up and down several times over a period of about 1 minute in each solution. A new transfer pipette should be used for the first wash solution but the same pipette can be used for subsequent transfers.
5. After the second wash solution, the oocytes are placed in the preferred culture medium equilibrated according to the manufacture recommendations at 37 °C for minimum 2 hours before fertilization.

STORAGE INSTRUCTIONS AND STABILITY
Store unopened vials refrigerated at 2-8 °C. Warm to room temperature (22-24 °C) or incubator temperature (37 °C), as appropriate, prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable in the unopened original vial until the expiration date shown on the label.
1. Remove desired volume of product using aseptic procedures. One vial contains adequate volume of medium for four thawing procedure.
2. Once removed, do not return any volume of product to the original vial. Discard remaining product.
3. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

Customer Service:
E-mail: customer.service@origio.com
Tel: +45 46 79 02 02
Fax: +45 46 79 03 02

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