Products:

Oocyte RS Single Lumen Needle

SL1603011	SL1603014
SL1703011	SL1703014
SL1753011	SL1753014
SL17530374B	SL1853011
SL1603011S	SL1603014S
SL1603014LS	SL1603064S
SL1603044S	SL1703011S
SL1703014S	SL1703064S
SL1703044S	SL1603511S
SL1703511S	

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Oocyte RS (Single lumen needle)**Oocyte RSS (Single lumen needle with stopcock)****Oocyte RD (Double lumen needle)**

The Oocyte RS, RSS and RD are retrieval needles for ultra-sound guided transvaginal aspiration and flushing of oocytes from ovarian follicles.

The products are for treatment of infertility. The products should only be used by professionals trained in IVF treatment.

WARNING:

- **Do not reuse. Reuse may cause contamination, patient infection, and failed procedure.**
- **Low resolution ultrasound equipment may compromise the visibility of the needle.**
- **When performing transabdominal ultrasound, ensure that the patient has a full bladder.**
- **To avoid infection only use sterile culture flushing media and sterile techniques. Infection may result in urinary tract infection, pelvic inflammatory disease, or uterine infection.**
- **When using our needle sets for the first time use as minimal force as possible until the user is familiar with the level of pressure they require in order to make a follicle puncture.**

Contraindications

The needles should not be used:

- In the presence of chronic cervical infection.
- In the presence of or after recent pelvic inflammatory disease.
- In the presence of a sexually transmitted disease.
- In the presence of or after a recent uterine perforation.
- In the presence of or after a recent caesarean section.
- In the presence of pregnancy.
- If the patient currently has an intrauterine device.
- For intra-fallopian procedures.

Cautions

- Read the complete Instruction For Use before use.
- For use of all tube length please refer to your SOP or clinical protocol for oocyte retrieval. Special precaution should be taken to follow your SOP or clinical protocol when using shorter tubing lengths.
- Ensure to use an appropriate needle guide with an inner diameter bigger than the outer diameter of the needle.
- Failure to properly follow the instructions, warnings and cautions may lead to serious surgical consequences or injury to the patient.
- Discard if product or packaging is damaged.
- Do not use the product if expiry date has been exceeded.

Quality control testing

Mouse Embryo Assay (MEA) tested.

Endotoxin (LAL) tested. (Ph. Eur., USP).

Note: The results of each batch are stated on a Certificate of Analysis, which is available on www.origio.com.

Storage instructions and stability

Store in original package at 0-30°C.

The product is sterilized using the gamma radiation processing technique and supplied sterile.

Radiation Indicator placed on the packaging turns red when exposed to radiation.

When stored as directed by the manufacturer the product is stable until the expiry date shown on the label.

Instruction for use

1. Place the patient in the appropriate position. Local anaesthetics may be administered as necessary.
2. Prepare the vulva and vagina for full aseptic technique and cover the area in sterile drapes.
3. Select the appropriate needle set. Selection of a single or double lumen needle should be based on a clinical need, training and experience.
4. Remove the needle from the packaging maintaining the sterility of the product.
5. The sterile needle set should be inspected for needle tip sharpness, kinked tubing & connections.
6. a) Attach the silicone bung to a sterile sampling tube.
b) Connect the luer lock on the vacuum tube of the needle set from the bung either directly to a luer compatible syringe, a disposable vacuum line or directly to pump.
c) Test the system to ensure all connections are sound and flow can be achieved.
d) Flush with appropriate media before use
i) Single lumen with Stopcock refer to Figure 1 'A'

WARNING

Failure to test the system may result in failure to obtain the required sample and/or tissue damage.

7. Attach the needle to the ultrasound transducer.
8. Place the ultrasound transducer and needle into the vaginal fornix to visualize the ovary and follicles. Identify the follicles to be aspirated.

WARNING

Failure to use ultrasound to locate the follicles may result in failure to obtain the required sample and/or tissue damage.

9. Under ultrasound guidance, advance the needle tip into the ovarian follicle. If the ovaries are mobile apply pressure suprapubically and enter the follicles with a stabbing movement. When the tip of the needle can be seen in the follicle, aspirate the content until empty using either the syringe or vacuum pump attached to the set. Gentle manipulation of the follicle may assist separation of oocytes from the follicle wall.
10. Using a vacuum pump, aspirate the follicle into the test tube.
 - a. i) Single lumen with Stopcock refer to Figure 1 'B'
11. If flushing is necessary, proceed with one of the following:
 - a) Single lumen: Detach the sampling tube from the bung. Attach a syringe with sterile media and flush down the aspiration tube. The aspiration tube protrudes 1 cm from the narrow end of the bung.
 - b) Single lumen with Stopcock: Attach a syringe with sterile media to the stopcock connection and ensure the stopcock valve is in the right position before flushing. See Figure 1 'C'
 - c) Double lumen: Attach a syringe with sterile media to the luer connection of the flushing tube and flush.
12. Repeat for all follicles within the ovary. To minimize trauma and possible bleeding keep the needle tip within the ovary between aspirations.
13. When all the follicles in one ovary have been aspirated, repeat the procedure with the other ovary.
14. Withdraw the needle from the patient, check the vagina for signs of bleeding and allow the patient to rest. If vaginal bleeding is experienced post collection, apply steady pressure with a sponge and forceps or use an absorbable suture if bleeding persists.
15. Discard the needle in accordance with local medical hazardous waste practices.

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