Enthält
Humanalbuminlösung (HAS)
Rekombinantes Humaninsulin
Gentamicinsulfat 10 µg/ml

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en - english

For slowing down the movement of the spermatozoa for ICSI.

This product is for ART treatment, whether the cause of infertility is male or female. The product should only be used by professionals trained in ART treatment.

Contains
Human albumin solution (HAS)
Recombinant human insulin
Gentamicin sulphate 10 µg/ml

Quality control testing
Sterility tested (Ph.Eur., USP)
Endotoxin tested ≤ 0.2 EU/ml (Ph.Eur., USP)
HSA analysis (Ph.Eur., USP)
Sperm motility tested
Sperm survival tested
Note: The results of each batch are stated on a Certificate of Analysis, which is available at www.origio.com.

Storage instructions and stability
The products are aseptically processed and supplied sterile.
Store in original container at 2-8 °C, protected from light.
Do not freeze.
The product is to be used within 7 days after opening.

When stored as directed by the manufacturer the product is stable until the expiry date shown on the label.

Precautions and warnings
Do not use the product if:
1. Product packaging appears damaged or if the seal is broken.
2. Expiry date has been exceeded.
3. The product becomes discoloured, cloudy, turbid, or shows any evidence of microbial contamination.

Caution: All blood products should be treated as potentially infectious. The source material used to manufacture this product was tested and found non-reactive for HBsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

Note: Due to the native physical characteristics of PVP, small variations in viscosity may occur between batches. The presence of human albumin solution (HAS) may occasionally cause a small amount of precipitation in the medium.

Note: Please note the need for traceability of this product. In addition, national legal requirements for this field may exist in your country.

Note: Only to be used in combination with other devices intended for the particular purpose.

Note: Dispose of the device in accordance with local regulations for disposal of medical devices.

Instructions for use
1. Remove PVP Clinical Grade and preferred holding medium from storage at 2-8°C and leave at room temperature for 10 minutes.
2. Depending on the number of oocytes for injection, pipette a corresponding number of 10 µl droplets of holding medium onto the bottom of the ICSI dish.
3. In the middle of the same dish place a 5-10 µl droplet of PVP Clinical Grade.
4. Cover with pre-equilibrated Liquid Paraffin and place the dish in a 5-6% CO₂ environment at 37°C for 30 minutes prior to use.
5. Introduce 2 µl of prepared and washed sperm to the droplet of PVP Clinical Grade. The PVP will reduce the motility of the sperm and facilitate the capture and loading of a single spermatozoon in the injection pipette.

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