English

Phosphate Buffered Saline (PBS)

For laboratory procedures only; other uses must be qualified by the end user.

	REF	
Product Description	Number	Unit Size
Phosphate Buffered Saline	ART-4011	500 mL

INTENDED USE

This product is intended for in vitro procedures involving the use of an isotonic solution with buffering capacity in the physiological range of vertebrate cells and tissues (pH 7.5 ± 0.2); (Dulbecco & Vogt, J Exp Med. 1954;99:167.).

PRODUCT DESCRIPTION

A phosphate buffered isotonic solution of the major inorganic salts required by cells.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter or cloudiness.

Single Use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled

expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

This solution is suggested for use in procedures that may benefit from the use of an isotonic, buffered solution of ions commonly found in physiological solutions.

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers at 20 to 25 °C. Warm to incubator (37 °C) temperature prior to use. Do not freeze or expose to temperatures greater than 40 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Call the SAGE™ SUPPORT LINE:
In the U.S.: (800) 243-2974
International: (203) 601-9818

EXPLANATION OF SYMBOLS

- REF** Catalog Number
- LOT** Batch Number
- Use By (year, month, day)
- Do Not Reuse
- Temperature Limitation
- Aseptic Technique Sterilization Membrane Filtered (SAL 10⁻³)
- ATTENTION:** See instructions for use.
- Authorized Representative in the European Community.
- Product conforms to the Medical Device Directive 93/42/EEC
- Manufacturer
- RX ONLY** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In Vitro Fertilization
a CooperSurgical Company

SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

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