English

**Quinn's Advantage™
Embryo Freeze Kit**

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Embryo Freeze Kit	ART-8014	5 x 12 mL

INTENDED USE
This product is intended for use in freezing pro-nuclear- and cleavage-stage embryos. This kit is designed to be used in conjunction with the SAGE™ Quinn's Advantage™ Thaw Kit (ART-8016) for warming and recovery of specimens.

PRODUCT DESCRIPTION
The components of this kit will allow for the efficient freezing of pronuclear- and cleavage-stage embryos. The components and recommended procedures are the preferred method for improved embryo survivability.

This product contains 10mg/L of genta-micin, an aminoglycoside antibiotic.

- MATERIALS PROVIDED IN THE EMBRYO FREEZE KIT**
- 2 x 12 mL vial of 1.5 M Propanediol, 0.1 M Sucrose Freezing Medium (REF #8001-12) with 12 mg/mL Human Serum Albumin
 - 2 x 12 mL vials of 1.5 M Propanediol Freezing Medium (REF #8003-12) with 12 mg/mL Human Serum Albumin
 - 1 x 12 mL vial of Freeze/Thaw Diluent Solution (REF #8013-12) with 12 mg/mL Human Serum Albumin

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PRECAUTIONS AND WARNINGS
Do not use medium that shows evidence of particulate matter, cloudiness.

This product contains albumin, a derivative of human blood.

The Quinn's Advantage™ 0.1 M Sucrose Freezing Medium (REF#8001-12) contains 12 mg/mL human serum albumin.

The Quinn's Advantage™ 1.5 M Propanediol Freezing Medium (REF#8003-12) contains 12 mg/mL human serum albumin.

The Quinn's Advantage™ Freeze/Thaw Diluent Solution (REF#8013-12) contains 12 mg/mL human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBSAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of

transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Single use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

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This product contains the antibiotic genta-micin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

QUALITY ASSURANCE
One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE FOR CRYOPRESERVATION OF EMBRYOS
Controlled hyperstimulation of women undergoing IVF or GIFT produces, on average, 10 to 12 mature oocytes for insemination. It is prudent to replace only a limited number of the resulting embryos, as multiple pregnancies can arise if too many embryos are replaced. Therefore, the majority of patients will have supernumerary embryos. These embryos can be cryopreserved and stored for later use, thus avoiding the necessity of the couple to undertake another stimulated cycle to recover more oocytes for IVF. The major cause of cell damage during cryopreservation is the formation of intracellular ice during freezing and thawing. By using cryopreservatives, controlling the rates of freezing and thawing, and carefully diluting the cryopreservative from the embryo after thawing, methods have been developed that allow 80% or more of frozen-thawed embryos to survive and be replaced into the reproductive tract of the woman who produced the oocytes or a genetically nonrelated recipient.

1. Prepare solutions containing 0.5 M and 1.0 M propanediol (PPD) by diluting the stock solution of 1.5 M PPD (REF #8003-12) with the Freeze/Thaw Diluent Solution (REF #8013-12):

- a. To prepare the 0.5 M PPD solution, add 0.3 mL of 1.5 M PPD (REF #8003-12) to 0.6 mL of diluent (REF #8013-12).
 - b. To prepare the 1.0 M PPD solution, add 0.6 mL of 1.5 M PPD (REF #8003-12) to 0.3 mL of diluent (REF #8013-12).
2. Embryos are pipetted at 37 °C into the 0.5 M PPD solution for 5 minutes, then the 1.0 M PPD solution for 5 minutes, and finally into 1 mL of the 1.5 M PPD solution (REF #8003-12) for 10 minutes. They are then transferred to 1 mL of 1.5 M PPD + 0.1 M Sucrose Freezing Medium (REF #8001-12) and pipetted into straws* containing this same solution. They are held at 37 °C in the 1.5 M PPD + 0.1 M Sucrose solution (REF #8001-12) for a total of 5 minutes before cooling is initiated.

3. As an alternative, the embryos can be placed directly into the 1.5 M PPD solution for 10 minutes before transfer to the 1.5 M PPD + 0.1 M Sucrose Freezing Medium.

*Vials, eg, 1.2 mL plastic cryovials, can also be used and have been found to be preferable by some laboratories.

It is important to make sure that the embryos are well mixed with the cryoprotectant solutions. This can be accomplished by pipetting the embryos up and down in the solution several times after adding them to the cryoprotectant solution. It is also recommended that the media be covered with Sterile Oil for Tissue Culture (REF #4008) during use to minimize evaporation of water and a subsequent change in osmolality of the solutions.

COOLING PROTOCOL
Embryos are usually frozen in straws, but 1.2 mL vials can be used as an alternative. The embryos are taken from a starting temperature of 37 °C to –6 °C at 2 °C/min. They are then seeded manually and held at –6 °C for a total

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of 10 to 15 minutes before being cooled at about 0.3 °C/min to around –35 °C. They are then transferred to a storage tank of liquid nitrogen.

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY
Store unopened containers at 2 °C to 8 °C upon receipt. Warm to incubator (37 °C) temperature prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label:

- A. Remove desired volume of product using aseptic procedures.
- B. Once removed, do not return any volume of product to the original container.
- C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

ART-8016 Quinn's Advantage™ Thaw Kit

ART-4008P Oil for Tissue Culture

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is trademark of CooperSurgical, Inc.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

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EXPLANATION OF SYMBOLS

- REF** Catalog Number
- LOT** Batch Number
-  Use By (year, month, day)
-  Do Not Reuse
-  Temperature Limitation
- STERILE A** Aseptic Technique Sterilization Membrane Filtered (SAL 10⁻³)
-  **ATTENTION:** See instructions for use.
- EC REP** Authorized Representative in the European Community.
- CE**  Product conforms to the Medical Device Directive 93/42/EEC
-  Manufacturer

RX ONLY U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In Vitro Fertilization
a CooperSurgical Company
SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

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Quinn's Advantage™ je ochranná známka společnosti CooperSurgical, Inc.