## Products:

### **Transem Spherical / EchoGen™ Spherical**

Product No.: TSET24SAS / TSET24SASEG

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The Transem EchoGen™ and Spherical catheters are for embryo transfer into the uterine cavity.

- The EchoGen™ range of catheters contains a metal band positioned at the distal end of the catheter tip enhancing visibility under ultrasound use.
- The Transem Spherical range of catheters has an atraumatic spherical tip at the distal end allowing reliable outer sheath position in difficult anatomical conditions.

The products are for treatment of infertility. The products should only be used by professionals trained in IVF treatment.

#### **Warning:**

- Do not reuse. Reuse may cause contamination, patient infection, and failed procedure.
- Low resolution ultrasound equipment may compromise the visibility of the catheter.
- When performing transabdominal ultrasound, ensure that the patient has a full bladder.
- To avoid infection only use sterile culture flushing media and sterile techniques. Infection may result in urinary tract infection, pelvic inflammatory disease, or uterine infection.
- To avoid infection ensure that the adjustable marker is still attached to the outer sheath after transfer.
- To avoid infection or loss never touch or adjust the metal band on the EchoGen range of catheters.

#### **Contraindications**

The catheters and stylet should not be used:

- In the presence of chronic cervical infection.
- In the presence of or after recent pelvic inflammatory disease.
- In the presence of a sexually transmitted disease.
- In the presence of or after a recent uterine perforation.
- In the presence of or after a recent caesarean section.
- In the presence of pregnancy.
- If the patient currently has an intrauterine device.
- For intra-fallopian procedures.

#### **Cautions**

- Read the complete Instruction For Use before use.
- Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.
- The transfer catheter must be used with the outer sheath.
- If additional forming of the catheter is needed, please note that it shall be done using thumb and forefinger. Before use, check that the outer sheath is not kinked and the inner catheter and obturator transitions smoothly through the outer sheath.
- Ensure that the adjustable cervical stopper and positioning band are still attached to the outer sheath after transfer.
- Discard if product or packaging is damaged.
- Do not use if expiry date has been exceeded.

#### **Quality control testing**

Mouse Embryo Assay (MEA) tested.

Endotoxin (LAL) tested. (Ph. Eur., USP).

Note: The results of each batch are stated on a Certificate Of Analysis, which is available on [www.origio.com](http://www.origio.com).

#### **Storage instructions and stability**

Store in original package at 0-30°C.

The product is sterilized using the gamma radiation processing technique and supplied sterile.

Radiation indicator dots placed on the packaging, turns red when exposed to radiation.

When stored as directed by the manufacturer the product is stable until the expiry date shown on the label.

### Instruction for use

1. Place patient in the appropriate position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.
2. When using transabdominal ultrasound follow normal procedures such that an optimal image of the uterine cavity is achieved.
3. The adjustable cervical stopper marker on the outer sheath is preset at 4cm from the distal end of the outer catheter tip and can be set to indicate the desired cervical length/depth. When moving this marker ensure that the orientation notch still corresponds with the direction of the pre-formed angle on the outer sheath.
4. The white positioning marker band set at 5cm from the distal end of the outer sheath and 1cm from the cervical stopper can be used to determine if the cervical stopper has moved during the embryo transfer procedure. If the cervical stopper does move check to ensure that the spherical tip of the outer sheath does not pass through the internal os. See Warning No. 2 below. The white positioning band can be adjusted by moving the cervical stopper as desired by the physician to determine the length / depth of the cervical canal.
5. **Two-stage transfer** (TSET24SASO / TSET24SASOEG): With the hubs of the obturator and outer sheath locked together, advance through the cervical canal to the inner os. If slight resistance is experienced, withdraw the catheter. Using thumb and forefinger curve the catheter with the obturator firmly locked in the outer sheath to complement the patient's anatomy and pass the catheter in this position to the internal os. Gently uncouple and remove the obturator leaving the outer sheath in place.

### WARNING

**- The hubs of the obturator and outer sheath must be locked together. Removing the obturator from the outer sheath prior to placement, may result in bleeding and damage to the endometrial tissue.**

6. Under sterile conditions, fill a 1 ml syringe with medium and eject any air bubbles. Remove the catheters tip protector and firmly attach the syringe to the luer compatible adaptor of the inner catheter. Inspect that the catheter is not kinked in the transition area between the catheter and the supporting catheter. Flush the catheter through with medium, remove syringe, and refill with medium. Reconnect the syringe and expel medium until the 0.1ml mark is reached, with the catheter tip pointing upwards to facilitate the removal of air bubbles.
7. Place the end of the catheter carefully into the medium droplet containing the embryos and inject a small amount of medium to break the boundary of surface tension that may appear.
8. Aspirate the embryos into the catheter and hand to the person performing embryo transfer.
9. **a) Single stage transfer** (TSET24SAS / TSET24SASEG: Inner catheter and outer sheath)  
With the hubs locked in position, advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity.
- b) Two-stage transfer** (TSET24SASO / TSET24SASOEG: Inner catheter and outer sheath)  
Advance the inner catheter through the outer sheath already positioned (continued from section 4. above) so that the inner

catheter passess through the external and internal os, into the mid-uterine cavity.

10.It may be necessary to twist the catheter as it negotiates the cervical canal. When using ultrasound, manipulation the probe to obtain the best image of the inner catheter.

### **WARNING**

- **To avoid bleeding and damage to the endometrial tissue:**

- 1. Never force the catheter against the digitally felt resistance while inside the uterine cavity.**
- 2. Never advance the outer sheath further than the internal os, and never enter the uterine cavity.**
- 3. Never pull back the inner catheter further than the last graduation marking.**
- 4. Never advance the tip of the catheter further than 1 cm from the fundus**
- 5. The supported inner catheter should not be used without the outer sheath.**

11.If slight resistance is experienced, withdraw the catheter and advance the outer sheath until only the tip of the inner catheter is visible. Align the most distal marking on the base of the inner catheter with the hub of the outer sheath should ensure a smooth radius tip protrudes. Whilst held in this position, pass the sheath and catheter to the internal os and then advance the inner catheter through the internal os, into the mid-uterine cavity. In the rare event where the cervical canal cannot be negotiated by this method stylets (Part No.: TSS1924) are available to assist insertion.

- Transem Stylet:

- a) Check/verify the perfect matching of the sheath to prevent the rigid rotation of the stylet within the malleable outer sheath.
- b) Pre-forming the stylet / sheath if necessary or required in order to facilitate the integration.
- c) Introduce the combined stylet / sheath in the cervical canal.
- d) Gently disengage and remove the stylet from the outer sheath and proceed with the embryo transfer procedure.

### **WARNING**

- **The warnings stated in section 10 are applicable for both catheters, stylets and obturators.**
- **To avoid perforation, special attention should be paid during the use of the Transem Stylet.**

12.Place tip of catheter approximately 1cm from the fundus. Twist the catheter to remove any mucus from the tip.

13.Expel the embryos.

14.Leave time for the uterine environment to equilibrate and slowly remove the catheter.

15.Hand catheter back to the embryologist to check that no embryos remain. Remove the speculum and allow the patient to rest.

16.Discard the catheter in accordance with local medical hazardous waste practices.

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