



AIR FILTRATION BY LIFEGLOBAL®

CodaAir® Systems

Keeping the lab air clean – the VOC removers



CodaAir Systems – mobile units for every lab

Unique VOC removing filter technology with a proven track record

- All products are clinically proven for safe and effective use in human ART laboratories
- Effective air filtration technology to purify VOCs (Volatile Organic Compounds) and pollutants
- High Performance™ 4-stage Filter® system has been proven to improve air quality and may increase expected outcome and overall results
- Coda systems have been used in human ART laboratories for more than two decades
- FDA cleared, CE certified and manufactured in the USA
- Patented technology; ISO 13485 certified

Product Code	Product Name	Quantity	Article Description
GCAU-115	Coda Aero® 500 - 115V	1 unit	Small, compact, versatile unit (300 - 600 sq.ft. lab space). Initial filters included
GCAU-230	Coda Aero 500 - 230V	1 unit	Small, compact, versatile unit (300 - 600 sq.ft. lab space). Initial filters included
GCAF-001	Coda Aero 400/500 Quarterly Filter	1 unit	Quarterly Filter (1 pk) - change filters every 3 months
GCAF-002	Coda Aero 400/500 Quarterly Filter	2 pack	Quarterly Filter (2 pk) - change filters every 3 months
A700-115	Coda Aero 700 - 115V	1 unit	Small, compact, versatile unit (500 - 700 sq.ft. lab space). Initial filters included
A700-230	Coda Aero 700 - 230V	1 unit	Small, compact, versatile unit (500 - 700 sq.ft. lab space). Initial filters included
GCAF-702	Coda Aero 700 Semi-Annual Filter	2 pack	Semi-Annual Filter (2 pk) - change filters every 6 months
A800-008	CodaAir 800 - 24V	1 unit	Small, compact, versatile unit (500 - 700 sq.ft. lab space). 24V system. Initial filters included
CAFK-800	CodaAir 800 - Semi-Annual Filter	2 pack	Semi-Annual Filter (2 pk) - change filters every 6 months
A900-009	CodaAir 900 - 24V	1 unit	Large capacity purification unit (900 - 1100 sq.ft. lab space). 24V system. Initial shirtbox filters included
CAFK-900	CodaAir 900 - Semi-Annual Filter	2 pack	Semi-Annual Filter (2 pk) - change filters every 6 months
A900-009	CodaAir 900 - 24V	1 unit	Large capacity purification unit (900 - 1100 sq.ft. lab space). 24V system. Add initial canister filter CCIK-900
CCIK-900	CodaAir 900 - Initial Filter Kit	1 kit	Initial Canister Filter (1 kit)
CCSK-900	CodaAir 900 - Semi-Annual Filter Kit	1 kit	Semi-Annual Filter (1 kit) - change filter kit 6 months after initial/annual kit
CCFK-900	CodaAir 900 - Annual Filter Kit	1 kit	Annual Filter (1 kit) - change filter kit 6 months after semi-annual kit
GELT-115	Coda ECO™ 1200 - 115V	1 unit	Large capacity purification unit (800 - 1800 sq.ft. lab space). 4-stage filters system. Add initial filter GEIK-1200
GELT-230	Coda ECO 1200 - 230	1 unit	Large capacity purification unit (800 - 1800 sq.ft. lab space). 4-stage filters system. Add initial filter GEIK-1200
GEIK-1200	Coda ECO 1200 - Initial Filter Kit	1 kit	Initial Filter (1 kit)
GESA-1200	Coda ECO 1200 - Semi-Annual Filter Kit	1 kit	Semi-Annual Filter (1 kit) - change filter kit 6 months after initial/annual kit
GEAK-1200	Coda® ECO™ 1200 - Annual Filter	1 kit	Annual Filter (1 kit) - change filter kit 6 months after semi-annual kit
GECO-115	Coda ECO 1500 - 115V	1 unit	Large capacity purification unit (1000 - 2400 sq.ft. lab space). 4-stage filters system. Add initial filter GEIK-1500
GECO-230	Coda ECO 1500 - 230V	1 unit	Large capacity purification unit (1000 - 2400 sq.ft. lab space). 4-stage filters system. Add initial filter GEIK-1500
GEIK-1500	Coda ECO 1500 - Initial Filter Kit	1 kit	Initial Filter (1 kit)
GESA-1500	Coda ECO 1500 - Semi-Annual Filter Kit	1 kit	Semi-Annual Filter (1 kit) - change filter kit 6 months after initial/annual kit
GEAK-1500	Coda ECO 1500 - Annual Filter Kit	1 kit	Annual Filter (1 kit) - change filter kit 6 months after semi-annual kit