









Keeping the box incubator air clean - the VOC remover



## Coda Incubator Units

## Unique VOC removing filter technology with a proven track record

- All products are clinically proven for safe and effective use in human ART laboratories
- Effective air filtration technology to purify VOCs (volatile organic compounds) and pollutants
- Contain a HEPA filter to remove particulate matter

- The Coda2 Unit has improved flow, which saves on energy, and is made of 100% recyclable materials
- FDA cleared, CE certified and manufactured in a clean room facility in the USA
- Patended technology; ISO 13485 certified







**Coda Unit** 

Coda SP

Coda 2 Unit

Product Code	Product Name	Quantity	Article Description
C2IU-010	Coda 2 Unit	1 unit	Place on top shelf of incubator
GCIU-010	Coda Unit	1 unit	Place on top shelf of incubator
GCSP-010	Coda SP	1 unit	Mount into the access hole of the incubator
GCIF-001	Coda Filters	1 unit	Change every 30 days

