

ICSI ORIGIO

Slimline Micropipettes

- Micropipettes since 1992
- No compromise on quality
- Customer-driven product development



"The design of the new ICSI Slimline pipet allows for an extremely gentle injection procedure, which practically eliminates oocyte damage."

Marlena Duke, M.Sc.,
Embryology Laboratory Director, RMA of New York

origio

a CooperSurgical Fertility Company

ICSI ORIGIO Slimline Injection Pipettes

Advantages



- Minimizes damage to the oocyte by creating a small entry point
- More parallel taper but still rigid when injecting
- Excellent control when puncturing the oolemma

Characteristics



- Short taper
- Long parallel wall section
- Small diameter

ICSI Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Bevel	Inner diameter	Other
ICSI SlimLine	MIC-SLM-XX	0, 20, 25, 30, 35, 45	0.5 mm	35-40°	4.3-4.9 µm	Parallel Taper
ICSI SlimLine	MIC-SLM-XX-B1.0	25, 30, 35	1 mm	35-40°	4.3-4.9 µm	Parallel Taper



CE
0086

Copyright ORIGIO a/s • Order No. 60010334 • Version 1 : April 20, 2017.
All information correct at time of print. Specifications are subject to change without notice or obligation on the part of the manufacturer.

ORIGIO a/s
Knardrupvej 2
2760 Måløv
Denmark
Tel: +45 46 79 02 00
www.origio.com



CooperSurgical Fertility Companies